
- DE** Bedienungsanleitung
- GB** Operating Instructions
- FR** Mode d'emploi
- NL** Handleiding
- IT** Istruzioni per l'uso
- ES** Instrucciones de uso
- PL** Instrukcja obsługi
- SK** Návod na použitie
- HR** Upute za rukovanje
- CZ** Návod k použití
- RO** Manual de utilizare
- BG** Указания за ползване

<b>DE</b>	<b>Bedienungsanleitung .....</b>	<b>4</b>
<b>GB</b>	<b>Operating Instructions .....</b>	<b>8</b>
<b>FR</b>	<b>Mode d'emploi .....</b>	<b>12</b>
<b>NL</b>	<b>Handleiding .....</b>	<b>16</b>
<b>IT</b>	<b>Istruzioni per l'uso .....</b>	<b>20</b>
<b>ES</b>	<b>Instrucciones de uso .....</b>	<b>24</b>
<b>PL</b>	<b>Instrukcja obsługi .....</b>	<b>28</b>
<b>SK</b>	<b>Návod na použitie .....</b>	<b>32</b>
<b>HR</b>	<b>Upute za rukovanje .....</b>	<b>36</b>
<b>CZ</b>	<b>Návod k použití.....</b>	<b>40</b>
<b>RO</b>	<b>Manual de utilizare .....</b>	<b>44</b>
<b>BG</b>	<b>Указания за ползване .....</b>	<b>48</b>

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<http://www.bresser.de/8851200>

### **RISK to your child.**

Aids with sharp edges and tips are sometimes used with this device. Please store the device and all of its accessories and aids out of the reach of children. There is a risk of INJURY.

This device contains electronic components that operate via a power source (batteries). Children should only use the device under adult supervision. Only use the device as described in the manual; otherwise, you run the risk of an ELECTRIC SHOCK.

Batteries should be kept out of children's reach. When inserting batteries please ensure the polarity is correct. Leaking or damaged batteries can cause injury if they come into contact with the skin. If you need to handle such batteries please wear suitable safety gloves.

Children should only use this device under supervision. Keep packaging materials (plastic bags, rubber bands, etc.) away from children. There is a risk of SUFFOCATION.

The chemicals and liquids provided should be kept out of reach of children. Do not drink the chemicals! Hands should be washed thoroughly under running water after use. In case of accidental contact with the eyes

or mouth rinse with water. Seek medical treatment for ailments arising from contact with the chemical substances and take the chemicals with you to the doctor.

### **FIRE-/ DANGER OF EXPLOSION!**

Do not expose the device to high temperatures. Use only the recommended batteries. Do not short-circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short-circuit, a fire or an explosion.

### **RISK of material damage.**

Never take the device apart. Please consult your dealer if there are any defects. The dealer will contact our service centre and send the device in for repair if needed.

### **TIPS on cleaning**

Remove the device from its energy source before cleaning (remove batteries)

Clean the exterior of device with a dry cloth. Do not use cleaning fluids so as to avoid causing damage to electronic components.

Clean the lens (objective) only with the cloth supplied or some other soft lint-free cloth (e.g. micro-fibre). Do not use excessive pressure - this may scratch the lens.

Protect the device from dust and moisture. The batteries should be removed from the unit if it has not been used for a long time.

## DISPOSAL

 Dispose of the packaging material/s as legally required. Consult the local authority on the matter if necessary.

 Do not dispose of electrical equipment in your ordinary refuse. The European guideline 2002/96/EU on Electronic and Electrical Equipment Waste and relevant laws applying to it require such used equipment to be separately collected and recycled in an environment-friendly manner. Empty batteries and accumulators must be disposed of separately. Information on disposing of all such equipment made after 01 June 2006 can be obtained from your local authority.

 In accordance with the regulations concerning batteries and rechargeable batteries, disposing of them in the normal household waste is explicitly forbidden. Please make sure to dispose of your used batteries as required by law – at a local collection point or in the retail market. Disposal in domestic waste violates the Battery Directive.

Batteries that contain toxins are marked with a sign and a chemical symbol.

- 1 battery contains cadmium
- 2 battery contains mercury
- 3 battery contains lead

## All Parts

1. Eyepiece
2. Focus knob
3. Objective turret
4. Stage
5. Mirror
6. Electronic light source
7. Base with battery compartment
8. Rubber base
9. Microscope arm
10. Microscope “Biotar”
11. Prepared Slides
12. Cover Plates
13. Blank Slides
14. Replacement bulb
15. Shrimp farming accessories
16. Microscope instruments
17. Micro-slicer
18. Hatchery

## Your Microscope

### What is a Microscope?

The microscope is made of two lens systems: the eyepiece and the objective. We’re presenting these systems as one lens each, so that it’s easier to understand. In reality, however, the eyepiece (1) and the objective in the turret (3) are made up of multiple lenses. The lower lens (objective) magnifies the

prepared specimen (e.g. 11), which leads to a magnified image of the specimen. The picture, which you can’t see, is magnified once more by the second lens (eyepiece, 1). Then, you can see the “microscope picture.”

### Assembly and Location

Before you start, choose an ideal location for using your microscope. On the one hand, it’s important there is enough light (normal observation). On the other hand, it’s important that you can quickly make the room dark when you want to (projector observation). Furthermore, I recommend that you place the microscope on a stable surface, since you won’t be able to get any satisfactory results on a shaky surface.

### Normal Observation

For normal observation, place the microscope in a bright location (near a window or desk lamp, for example). Remove the microscope from the package and tilt the microscope arm (9) so that it’s in a comfortable viewing position for you.

The focus knob (2) is turned to the upper stop and the objective turret (3) is set to the lowest magnification.

Now, take a look through the eyepiece and adjust the mirror (5) so that you get a uniformly bright circle of light. Or, you can use the lamp. You’ll find further tips about the lamp

in the next chapter. Now, place a prepared slide (e.g. 11) under the clips on the stage (4), directly under the objective. When you take a look through the eyepiece, you can see the magnified specimen. At this point, it might still be a slightly fuzzy picture. The sharpness is adjusted by slowly turning the focus knob (2). You can now select a higher magnification by turning the objective turret and selecting a different objective.

In doing so, note that the sharpness of the picture must be adjusted again for the higher magnification. Also, the higher the magnification, the more light you will need for good illumination of the picture.

**TIP:** Avoid positioning the mirror so that it is exposed to direct sunlight, since it can cause a glare, which won’t allow you to get a clear picture.

### Electric Illumination

To use the electric light source (6), you need two 1.5 V batteries, which are installed in the base of the microscope (7). Remove the rubber base (8) on the microscope and insert the batteries according to the displayed +/- information. Reattach the rubber base to the microscope.

You can turn on the light by rotating it in the direction of the stage (4). Take a look through the eyepiece and adjust the light so that you get optimal picture brightness. Now you can observe in the same way as described in point 1.3.

The bulbs for the electric light can be replaced. One replacement bulb (14) is included. When you use another bulb, make sure to pay attention to the maximum wattage, which is displayed on the lamp body.

**TIP:** The higher the magnification used, the more light is required for a good illumination of the picture. Therefore, always start your experiments with one of the lower magnifications.

### Warranty & Service

The regular guarantee period is 2 years and begins on the day of purchase. To benefit from an extended voluntary guarantee period as stated on the gift box, registration on our website is required.

You can consult the full guarantee terms as well as information on extending the guarantee period and details of our services at [www.bresser.de/warranty\\_terms](http://www.bresser.de/warranty_terms).

Would you like detailed instructions for this product in a particular language? Then visit our website via the link below (QR code) for available versions.

Alternatively you can also send an email to [manuals@bresser.de](mailto:manuals@bresser.de) or leave a message on +49 (0) 28 72 – 80 74-220\*. Please always state your name, precise address, a valid phone number and email address, as well as the article number and name.

\*Number charged at local rates in Germany (the amount you will be charged per phone call will depend on the tariff of your phone provider); calls from abroad will involve higher costs.

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**Bresser GmbH**